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
What is DNA Extraction?

• (DNA) extraction is the process by which DNA is separated from proteins, membranes, and other cellular material in the cells

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#### Step 1: Lysis.

a. (Break the cells open to expose DNA)

b. Remove membrane lipids by adding detergent

#### Step 2: Precipitation

a. (Precipitate DNA with an alcohol — usually ethanol or isopropanol. Since DNA is insoluble in these alcohols,

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
works every time

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#### In this activity.

- 1-students crush strawberries and use detergent and salt to break the cell wall so that DNA will be released into the solution.
- 2- Alcohol is then poured slowly over the top of the solution and DNA will rise to the surface.
- 3- Students can spool the DNA with a stirring rod. Because strawberries have multiple copies of DNA (octoploid).

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## Procedure

- 1- Obtain one fresh or one frozen and thawed strawberry. If you are using a fresh strawberry, remove the green sepals (tops) from the berry.
 

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٦. اهرسي الفراولة لمدة دقيقة إضافية.
- 7-Place a funnel into a 50-mL centrifuge tube. Fold the cheesecloth in half along the longer side and place it in the funnel to create a filter. The cheesecloth will overlap the edge of the funnel.
 

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10-Hold the 15-mL tube at an angle. Using a transfer pipet, carefully add 5 mL of cold 95% ethanol by running it down the inside of the tube. Add the 95% ethanol until the total volume is 7 mL (use the lines on the side of the tube to help you measure). You should have two distinct layers.

Caution: Do not mix the strawberry extract and the ethanol!

11-Watch closely as translucent strands of DNA begin to clump together where the ethanol layer meets the strawberry extract layer. Tiny bubbles in the ethanol layer will appear where the DNA precipitates.

12- Slowly and carefully rotate the wooden stick in the ethanol directly above the extract layer to wind (or "spool") the DNA. Remove the wooden stick from the tube and observe the DNA .

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